

## Erythropoietin (EPO)

### Luminescence Cell-Based Potency Assay for Erythropoietin

Cell-based potency determination using BAF3 cells with overexpressed huEpoR receptors. ICH Q2(R2)-qualified methodology measuring ATP production after receptor activation. Suitable as 3R-compliant alternative for batch release and stability evaluation under GMP conditions.

#### WHO Standard-Based Methodology

- **Reference Standard:** WHO Standard Erythropoietin (NIBSC Code: 11/170)
- **Alternative to:** Ph. Eur. monograph 1316 Erythropoietin concentrated solution (ECS)

#### Ahead of Regulatory Trends



ICH Q2(R2)-qualified methodology supports your workflows for emerging acceptance of in vitro alternatives.

#### Screen Smarter, Save Faster



Eliminate upfront animal testing costs for candidate screening and early validation, saving upfront on in vivo models.

#### WHO Standard Foundation



7.16% mean GCV precision using NIBSC 11/170 reference standard accelerates your path to in vivo replacement or biosimilar submissions.

#### Technical Specifications

- Luminescence detection (RLU)
- 10 dose steps, triplicates per step
- 1 test sample + assay control
- 4-PL fit for potency determination
- USP <1032> equivalence margins



#### Applications

- Biosimilar characterization
- Batch release and stability testing
- Method qualification per ICH Q2(R2)
- Regulatory submission support
- Clinical trial material testing

# Potency Assay Example

## Assay Setup and Evaluation

- WHO Standard Erythropoietin (NIBSC Code: 11/170) was used with 50% and 150% potency level as sample
- A 4-PL fit was used for potency determination
- Titration was performed over 10 dose steps in triplicate
- Assays “passed” according to Equivalence Margins (USP <1032>)
- Relative Light Units (RLU) are depicted on graph as response
- Erythropoietin concentration depicted on graph as relative dose (log base 2)




# Method Qualification


## Result Summary

- Five potency level investigated (n = 6 each)
- Assessment of accuracy, linearity, precision (overall and intermediate), and repeatability

### Accuracy



### Linearity




## Precision (Overall)

RP Level	Geometric Mean	%GCV
0.50000	0.48453	10.63
0.75000	0.71390	7.05
1.00000	0.95539	5.08
1.25000	1.21361	7.38
1.50000	1.55188	5.69
	Mean	7.16

## Precision (Intermediate)

Mean	0.96624	GM	0.96431
SD	0.068	GSD	0.069
CV%	7.04	GCV%	7.13

## Repeatability (n = 6, Analyst 1)



**Let's Discover Together.**

Position your EPO program ahead of the curve.  
Speak with an expert.

