

# Bioanalytical Service

## Cell based Potency Assays

### G-CSF potency assay

- ➲ Determination of biological activity of G-CSF (Filgrastim)
- ➲ Validated method acc. to the European Pharmacopoeia (Monograph no. 2206)
- ➲ GMP-certified execution
- ➲ NFS 60 cell line
- ➲ Measurement of tetrazolium conversion

### G- CSF potency assay example

- ➲ Calculated 80 % -activity
- ➲ Standard 100 % (640 pg/ml)
- ➲ Titration factor 1.5
- ➲ Measurement: 77.8 % activity
- ➲ Confidence interval: 7.0 %

